

# Quick Reference Guide

## HORUS SCOPE 1 on AVIZIA CART

### Step 1

- ❖ Ensure the device is adequately charged (solid blue light on power indicator/handle when the device is powered on).



### Step 2

- ❖ Ensure camera/device is connected to cart: USB and Composite cable
  - Use connections in drawer of Avizia cart to connect camera to cart:
    - Mini AV cable from USB charger in drawer to bottom of camera.
    - Insert composite cable to camera and attach to port in inside drawer (yellow end on cable).



### Step 3

- ❖ Attach appropriate viewing lens as follows:
  - Ensure the device is powered off.
  - Find and match the lens alignment marks (small black circles).
  - Insert the lens into the Horus Scope connector interface.
  - Rotate the lens attachment gently clockwise until it stops.
  - Ensure the white position bars are in the appropriate horizontal position.



### Step 4

Press the power button to begin using the camera.

- ❖ Press "Aux" button on front of Avizia cart (below monitor) to share image(s)
- ❖ Press back button on remote to get out of preview screen (showing wireless and/or PC option): will then see image on screen/monitor



## Step 5

- ❖ Adjust brightness/focus as needed:
  - Brightness controls on the front (under “OK” button) adjust amount of LED light emitted from the device.
  - Brightness controls also act as zoom in/out buttons for still images.
  - Focus wheel (found on back) adjusts image quality.



## Step 6

Press “OK” button to freeze an image.



## Step 7

To stop sharing image, press “Aux” to return to video call with other participants.

## Step 8

- ❖ After call is ended, Power off the device.
- ❖ **Disinfect** the device handle with an alcohol wipe.
- ❖ Do not sterilize the device.
- ❖ Clean the lens with **dust free** cleaning cloth.



## Step 9

- ❖ Recharge device by connecting mini AV to USB in cart drawer (power indicator light will turn off indicating device is fully charged while connected to charger).
- ❖ Flashing blue light (on power indicator/handle) means battery is less than 25% full when unit is on.
- ❖ Orange light means battery is less than 100% full when connected to charger.



## Tips and Notes

---

- Horus Scope will not display an image if there is no lens attached.
- Do not touch the lens.
- Images/videos are not stored/recorded during virtual care appointments.
- To maintain charge, keep camera plugged in to USB port.
- If issues with camera, email [virtualcaresupport@nlchi.nl.ca](mailto:virtualcaresupport@nlchi.nl.ca) or call 1-888-590-1111